

# Lab meeting 7/28/09





What did the hippocampus say during its retirement speech?


**Thanks for the memories...**


- Development of a longer yet stable CAG repeat in mouse model:
  - Modeling a more toxic human pathology
    - Earlier onset of symptoms

### Q127

- Human SCA-2 driven by the PCP2 promoter
  - Transgene should be confined to the Purkinje Cell







Cummings et al, 2001; 10 (14): 1511-18 Human Molecular Genetics

# Methods: brain preparation

- Brain perfused with ice-cold PBS
- Overnight in 4% PFA at 4c
- Cryoprotection: 20% $\rightarrow$ 30% sucrose
- Embedded within OCT medium
  - Brain positioned for sagittal slices of the vermis
- Kept at -20c until time of slicing


# Methods: slice

- Cryostat set to -20c
- 10 $\mu$ M sections – slide mounted
- Stored at -80c until time of staining


# Methods: Stain

- 3 x 10min 1 x PBS wash
- Permeabilize (0.3% Triton) and block (10% NDS) in 1 x PBS @ room temp, 2hrs
- Primary antibody (714- 15 $\mu$ g/mL); Calbindin ab (1:300) in 0.1% Triton – 5%NDS in 1 x PBS overnight at 4c
- 3 x 15min 0.1% triton and 1% NDS in 1 x PBS
- Secondary ab: Cy3 (red; anti-mouse), Cy2 (green; anti-rabbit)
- 3 x 10min wash: 0.1% triton and 1% NDS in 1 x PBS (dapi 1:6k on final wash)
- 3 x 5min 1 x PBS
- Cover with slip with Prolong-Gold → dry overnight

All camera settings are held  
constant across slides



Sca2 Ko -/-




Sca2 Ko WT

Sca2Ko: 23 -/-  
DOB: 1/8/09  
10  $\mu$ M; 10x

Calbindin: red  
Sca2 ab (714): green  
Dapi: blue


Control Sca2 KO



GL PL ML



**Sca2Ko: 23 -/-**  
DOB: 11/3/08  
10  $\mu$ M; 10x



Calbindin: red  
Sca2 ab (714): green  
Dapi: blue






Sca2Ko: 24 -/-  
DOB: 11/3/08  
10  $\mu$ M; 10x

Calbindin: red  
Sca2 ab (714): green  
Dapi: blue






Sca2Ko: 24 -/-  
DOB: 11/3/08  
10  $\mu$ M; 10x



Calbindin: red


40-2: 6 months



Q127 40-2: 105 (wt)  
DOB: 12/8/08  
10  $\mu$ M; 10x

Calbindin: red  
Sca2 ab (714): green  
Dapi: blue








Q127 40-2: 100 (+)

DOB: 12/8/08

10  $\mu$ M; 10x



Calbindin: red  
Sca2 ab (714): green  
Dapi: blue





**Q127 40-2: 100 (+)**  
DOB: 12/8/08  
10  $\mu$ M; 10x


Calbindin: red  
Sca2 ab (714): green  
Dapi: blue




**Q127 40-2: 100 (+)**  
DOB: 12/8/08  
10  $\mu$ M; 10x

Calbindin: red  
Sca2 ab (714): green  
Dapi: blue

Arrows indicate heterotopically  
Located PCs



Arrows indicate heterotopically  
Located PCs



## 40-2: 8 weeks

- The following brain slices were not processed with the slides previously shown
- Calbindin stain only
- Coronol sectioning


**Q127 40-2: 139 (+)**

DOB: 5/8/09

10  $\mu$ M; 10x

Stained 7/21/09

Calbindin: red




**Q127 40-2: 135 (wt)**

DOB: 5/8/09

10  $\mu$ M; 10x

Stained 7/21/09

Calbindin: red






# Future

- IHC profile across age & genotype
- WB and Q PCR of same
- Test
- Breed Q127 animals with existing Park KO animals
- Ca++ functionality

# Extra slides

# Q127: 40-2: rotarod





**Q127 40-2: 105 (wt)**  
DOB: 12/8/08  
10  $\mu$ M; 10x

Calbindin: red  
Sca2 ab (714): green  
Dapi: blue

